



FLUORIDE in URINE

8308

F

MW: 19.00

CAS: 16984-48-8

RTECS: LM6290000

METHOD: 8308, Issue 3

EVALUATION: FULL

Issue 1: 15 February 1984

Issue 3: 14 March 2016

BIOLOGICAL INDICATOR OF: exposure to inorganic fluorides [1,2]

SYNONYMS: None

SAMPLING	MEASUREMENT
<p>SPECIMEN: urine, pre- and post-shift</p> <p>VOLUME: 50 mL in chemically clean polyethylene bottles</p> <p>PRESERVATIVE: 0.2 g EDTA added to bottles before collection</p> <p>SHIPMENT: in insulated containers using bagged refrigerant</p> <p>SAMPLE STABILITY: 2 weeks @ 4 °C, longer if frozen</p> <p>CONTROLS: collect 3 sets of specimens from unexposed workers (pre- and post-shift)</p>	<p>TECHNIQUE: ION SELECTIVE ELECTRODE (ISE)</p> <p>ANALYTE: fluoride ion (F⁻)</p> <p>DILUTION: mix equal volumes of urine with TISAB</p> <p>CALIBRATION: solutions of sodium fluoride in water</p> <p>QUALITY CONTROL: spiked urine pools; correct for creatinine content</p> <p>RANGE: 1 to 100 mg/L urine</p> <p>ESTIMATED LOD: 0.1 mg/L urine</p> <p>RECOVERY: 0.95 [3]</p> <p>PRECISION (\bar{S}_r): 0.04</p>
ACCURACY	
<p>RANGE STUDIED: not studied</p> <p>BIAS: not determined</p> <p>OVERALL PRECISION (\bar{S}_{rT}): not determined</p> <p>ACCURACY: ± 23.6%</p>	

APPLICABILITY: Any fluorine-containing substances that can be metabolized to fluoride (F⁻) can be monitored using this procedure. Inorganic compounds of fluoride can be absorbed by the body resulting in the excretion of fluoride ions as sodium fluoride. Dietary and domestic water sources of fluoride must be considered, as well as dental treatments.

INTERFERENCES: Hydroxide, the only positive interference, is eliminated by use of the buffer. Negative interferences from complexation of fluoride by cations, such as calcium, are minimized by EDTA preservative and the high ionic strength buffer.

OTHER METHODS: This method is P&CAM 114 [4] in a revised format. Other methods that have been used are those described in the NIOSH criteria documents on inorganic fluorides [1] and hydrogen fluoride [2].

REAGENTS:

1. Distilled or deionized water.
2. Sodium citrate tribasic dihydrate ($\text{Na}_3\text{C}_6\text{H}_3\text{O}_7 \cdot 2\text{H}_2\text{O}$), ACS reagent grade or better.
3. Ethylenediaminetetraacetic acid (EDTA) disodium salt, ACS reagent grade or better.
4. Acetic acid, glacial, ACS reagent grade or better.
5. Sodium chloride, ACS reagent grade or better.
6. Sodium hydroxide, 5 M. Dissolve 20 g NaOH in distilled water; dilute to 100 mL.
7. Sodium fluoride, ACS reagent grade or better.
8. Calibration stock solution, 100 $\mu\text{g F}^-/\text{mL}$. Dissolve 0.2211 g dry sodium fluoride in distilled water. Make 1000 mL solution.
9. Total ionic strength activity buffer (TISAB), pH 5. Add 57 mL glacial acetic acid, 58 g sodium chloride, and 0.30 g sodium citrate to a 1-L beaker containing 500 mL distilled water. Stir to dissolve. Place beaker in water bath for cooling. Slowly add 5 M sodium hydroxide until the pH is between 5.0 and 5.5. Cool to room temperature; dilute to 1 L with distilled water.
10. Fluoride in Urine Certified Reference Materials (CRMs).

EQUIPMENT:

1. Polyethylene bottles, 125-mL, wide-mouth.
2. Fluoride ion specific electrode (ISE), with reference electrode.
3. pH/millivolt meter, reading to ± 0.5 mV.
4. Stirrer, magnetic.
5. Stirring bars, PTFE-coated.
6. Beakers, plastic, 50-mL.
7. pH electrode.
8. Pipets, appropriate sizes for standards.
9. Volumetric flasks for standards.
10. Water bath.
11. Tissues, low-lint lab wipers.

SPECIAL PRECAUTIONS: Wear gloves, lab coat, and safety glasses while handling all chemicals and human urine products. Disposable plastic, glass, and paper (pipet tips, gloves, etc.) that contact urine should be placed in a biohazard container. Standard precautions should always be used when handling bodily fluids and/or extracts of bodily fluids [5]. Handle urine samples and urine extracts using proper gloves. Glacial acetic acid is flammable and corrosive and should be handled in a fume hood.

SAMPLING:

1. Collect pre- and post-shift spot urine samples in polyethylene bottles containing 0.2 g EDTA.
2. Ship samples in insulated container at about 4 °C using bagged refrigerant.

SAMPLE PREPARATION:

3. Perform a creatinine determination on an aliquot of the urine (e.g., [6]).

CALIBRATION AND QUALITY CONTROL:

4. Prepare at least five working standards in the range 0.1 to 100 $\mu\text{g F}^-/\text{mL}$ by appropriate dilutions of the calibration stock solution with distilled water.

5. Analyze a set of working standards together with the samples and blanks (steps 9 through 12) starting with the lowest concentration.
NOTE: Working standards, samples, and blanks must be analyzed under the same conditions, including temperature, for accurate results.
6. Prepare a semi-log calibration curve plotting millivolts on the linear scale (y-axis) and fluoride concentration, $\mu\text{g/mL}$, on the log scale (x-axis).
7. Maintain standardization by running a standard with every 10 specimens.
8. Run a spiked urine control specimen with every 10 specimens to maintain quality assurance.
NOTE: Urine used for spiked controls must be analyzed before use to determine background fluoride concentration.

MEASUREMENT:

9. Add 10 mL well-mixed urine and 10 mL TISAB to a 50-mL plastic beaker.
10. Place a small stirring bar into beaker and mix continuously on a magnetic stirrer at room temperature.
11. Immerse electrodes. Allow sample to mix for 2 to 3 min and then record millivolt reading.
12. Rinse electrodes and stirring bar thoroughly with distilled water and wipe dry with tissue before next sample analysis.

CALCULATIONS:

13. Convert the millivolt readings to fluoride concentration using the calibration curve from step 6.
14. Express fluoride concentration as mg F/g urinary creatinine.

GUIDES TO INTERPRETATION:

Urine concentrations of fluorides in normal non-occupationally exposed workers have been reported to range from 0.2 to 3.2 mg/L depending on dietary intake [7]. Pre-shift levels of less than 4 mg/g creatinine and post-shift levels of less than 7 mg/g creatinine appears to protect workers against bony fluorosis [8]. NIOSH has recommended that post-shift urine specimens should not exceed 7 mg/L (corrected to a specific gravity of 1.024) and pre-shift specimens should not exceed 4 mg/L (corrected to a specific gravity of 1.024) [1,2].

The Biological Exposure Indices (BEI) for fluoride (as of the date of this method's publication) are 2 mg/L prior to shift and 3 mg/L at end of shift [9]. This BEI changed in 2011.

EVALUATION OF METHOD

No formal method evaluation has been reported; however, Tusl [3] reported recoveries of added fluoride from 94 to 100%. Precision based on analysis of 25 specimens in triplicate is estimated to be better than $\bar{S}_r = 0.04$. This method employs standard methodology that has been shown to provide adequate performance data for decades. Additional evaluation data may be found in, but is not limited to, the following references [10-13].

REFERENCES:

- [1] NIOSH [1976]. NIOSH criteria for a recommended standard: occupational exposure to inorganic fluorides. Cincinnati, OH: U.S. Department of Health, Education, and Welfare, Center for Disease Control, National Institute for Occupational Safety and Health, DHEW (NIOSH) Publication No. 76-103.
- [2] NIOSH [1976]. NIOSH criteria for a recommended standard: occupational exposure to hydrogen fluoride. Cincinnati, OH: U.S. Department of Health, Education, and Welfare, Center for Disease Control, National Institute for Occupational Safety and Health, DHEW (NIOSH) Publication No. 76-143.

-
- [3] Tusi J [1970]. Direct determination of fluoride in human urine using fluoride electrode. *Clin Chim Acta* 27:216-218.
- [4] NIOSH [1977]. Fluoride in urine: P&CAM 114. In: Taylor DG, ed. NIOSH manual of analytical methods. 2nd ed., Vol 1. Cincinnati, OH: U.S. Department of Health, Education, and Welfare, Center for Disease Control, National Institute for Occupational Safety and Health, DHEW (NIOSH) Publication No. 77-157A.
- [5] CDC [2007]. 2007 Guidelines for isolation precautions: preventing transmission of infectious agents in healthcare settings [<http://www.cdc.gov/hicpac/2007IP/2007isolationprecautions.html>]. Date accessed: March 2016.
- [6] Tietz NW [1976]. *Fundamentals of clinical chemistry*. 2nd ed. Philadelphia, PA: W. B. Saunders Co., pp. 994-999.
- [7] Baselt RC [1980]. *Biological monitoring methods for industrial chemicals*. Davis, CA: Biomedical Publications, pp. 140-143.
- [8] Lauwreys RR [1983]. *Industrial chemical exposure: guidelines for biological monitoring*. Davis, CA: Biomedical Publications, pp. 26-27, 134.
- [9] ACGIH [2014]. *TLVs and BEIs based on the documentation of the Threshold Limit Values for chemical substances and physical agents and Biological Exposure Indices*. Cincinnati, Ohio: American Conference of Governmental Industrial Hygienists [www.acgih.org]. Date accessed: May 2014.
- [10] Singer L, Armstrong WD [1973]. Determination of fluoride in ultrafiltrates of sera. *Biochem Med* 8:415-422.
- [11] Chiba K, Tsunoda K, Haraguchi H, Fuwa K [1980]. Determination of fluoride in urine and blood serum by aluminum monofluoride molecular absorption spectrometry and with a fluoride ion selective electrode. *Anal Chem* 52:1582-1585.
- [12] Milde D, Nováková K, Čermáková I [2004]. Fluoride determination in urine with fluoride ion selective electrode: within laboratory method and sample storage optimization. *Acta Univ Palacki Olomu, Fac Rerum Nat, Chem* 43:104-109.
- [13] Singh B, Gaur S, Garg VK [2007]. Fluoride in drinking water and human urine in Southern Haryana, India. *J Hazard Mater* 144:147-151.

METHOD WRITTEN BY:

William P. Tolos, NIOSH.

Disclaimer: Mention of any company or product does not constitute endorsement by NIOSH. In addition, citations to websites external to NIOSH do not constitute NIOSH endorsement of the sponsoring organizations or their programs or products. Furthermore, NIOSH is not responsible for the content of these websites. All web addresses referenced in this document were accessible as of the publication date.