

# Jasco Digital Polarimeter Instructions

Printed: September 24, 2021

These instructions will show the user how to collect Optical rotation data on their samples. Currently we have two sample cells, one has a 50mm path length and the other has 100mm path length. Either of these cells may be checked out from me. The instrument is set to a frequency of 589nm. Other values are available but filters and lamps will need to be purchased.

In this instruction set, I will be using the following conventions:

- Italicized*    Typed keyboard entry in appropriate window
- Boldface**    Click on the item with the left mouse button  
                  Item may be a menu entry or button on dialog window

## OPERATING INSTRUCTIONS

### LOG ON TO INSTRUMENT

Move mouse to activate screen

LOGIN: *Username*                      Log User into the Instrument Facility system

PASSWORD: *Password*

**OKAY**

### START HARDWARE AND SOFTWARE

Push the On/Off button              Lower left corner of Console

**SpectraManager**                      Single Click              Program starts

**Standard Measurement**              Double Click              Start data collection module

If needed

**Polarimeter Analysis**              Double Click              Start analysis module

# Jasco Digital Polarimeter Instructions

## COLLECT DATA

<b>Measure</b>	Select Measure in the Menu
<b>Parameters</b>	Select parameters from Dropdown Menu
Change values as needed	
Source	W1 Do not change
D.I.T	Length of Sample time
Cycle Time	Numbers of reads to take
Cycle Interval	Time between reads if Cycle Time > 1
<b>B ICON</b>	Collect Baseline sample – Sample cell should hold only solvent
<b>Measure</b>	Start measurement
<b>S ICON</b>	Collect Sample Data
<b>Sample Name</b>	Enter Sample name
<b>Comment</b>	Enter Comment
<b>Measure</b>	Start measurement
Start Button	Same as S-ICON but does not ask for input

## SAVE DATA

<b>File</b>	Select File in the Menu
<b>Save(as)</b>	Save data in proprietary format (Read only by this software) Data may be reviewed using Polarimeter Analysis program
<b>Export</b>	Save data as ASCII text file (Read by other programs such as EXCEL)

## SHUTDOWN and EXIT

<b>File</b>	Select File in the Menu
<b>Exit</b>	Select Exit from Dropdown Menu
<b>Program</b>	Select Program in the Menu
<b>Exit</b>	Select Exit from Dropdown Menu
Push the On/Off button	Lower left corner of Console
<b>Log Out</b>	Log off instrument

# **Jasco Digital Polarimeter Instructions**

## **SIGN UP POLICY for the Olis Globalworks Cary 14 Spectrophotometer**

This instrument is available on a walk in basis. There is currently no signup policy.

## **DISK POLICY for the Jasco Digital Polarimeter**

Disk Storage: 35.0 GBytes Available Storage

This data is mirrored on the NMRNET Server every 15 minutes

Data will remain on this disk for 90 days after which it will be deleted

Server Storage: 12.5 GBytes Total Storage

The data will stay on the server for 90 days after which it will be either deleted or sent to an archive location at the group advisor's decision.

Your data may be retrieved from the Facility Server using the FTP program on your computer. The parameters are:

Hostname: nmrnet ([nmrnet.usc.edu](http://nmrnet.usc.edu) for off campus access)

Username: 'Group name'

Password: 'Group password'

Data: nmrdata\*username*\POLARIM

Make sure to do a binary transfer.