

Using CPX3 to extract data

- (1) Start up cpx3
- (2) Switch to the desired data mode (raw, edited, clean, average) using the 'display' pull down menu.
- (3) Use the 'File' pulldown menu to chose 'export data'
- (4) Set the time range and select variables from the list
- (5) Click on 'output'
- (6) Change to the convert/process data tab. Eventually will get a screen about exporting.
- (7) Check the box the 'file' box.
- (8) Click on disk icon to set file name and directory. Click save. This will put file name and path in text window
- (9) Choose file option settings (they are set to common defaults so you may not need to)
 - a. Each option choice opens a pulldown window or checkbox for that choice
 - b. Make all file option choices
- (10) Once all file option choices have been made click 'OK' and the file will be generated.

TIP1: You can see what the file will look like if you skip steps 7&8. In this case, once you click 'OK' in step 10, a window will pop up with the data. We recommend using file viewing option only for shorter time series (e.g., 1 week).

TIP2: Once you are happy with how the file looks in the viewing window then extend the time range (step 4) and do steps 7&8 to generate the final file.