



I seek the future.

**MiSeq.** Next-generation sequencing for all you seek.

You want amazing accuracy and performance on your benchtop. Illumina delivers—yet again. The MiSeq Personal Sequencer is the only fully integrated, truly end-to-end benchtop solution around. It's just one more example of why Illumina solutions generate a remarkable 90% of all the world's sequencing data. Discover what's possible.

[www.illumina.com/iseek](http://www.illumina.com/iseek)

**illumina**<sup>®</sup>



I AM AN FFPE SPECIMEN. SCIENTISTS  
HOPE TO MAKE SIGNIFICANT

# PROGRESS

IN UNDERSTANDING MY DISEASE STATE  
THROUGH ADVANCED BIOINFORMATICS.  
IT'S MY TURN TO CONTRIBUTE TO THE CURE.  
THAT'S WHY I'M ON THE WAY TO  
EXPRESSION ANALYSIS.

EVERY SPECIMEN HAS A STORY. LET YOURS BE HEARD AT [GENOMICKNOWHOW.COM](http://GENOMICKNOWHOW.COM)



GENOMIC KNOW-HOW®



PGM™ for genes. Proton™ for genomes.  
**Sequencing for all.**

Powered by fast, simple, scalable semiconductor chips, the Ion PGM™ Sequencer introduced an entirely new approach to sequencing, making it dramatically faster and more accessible.

The new Ion Proton™ Sequencer will go even further. With chip densities up to 1,000 times greater than the Ion PGM™ Sequencer, the Ion Proton™ Sequencer will put whole-genome sequencing within reach of every lab.

Find fast, affordable benchtop sequencing at  
[lifetechnologies.com/ionsequencing](http://lifetechnologies.com/ionsequencing)



# Directional RNA-Seq Libraries in 4 Hours.

From just 500 pg RNA.

Get RNA-Seq libraries on the sequencer the same day, with the new ScriptSeq™ v2 Kits. Start with just 500 pg of rRNA-depleted or poly(A)-enriched RNA. No cDNA fragmentation, adaptor ligation, or gel size-selection.

Go ahead. Make *your* day.

[www.epicentre.com/scriptseqv2](http://www.epicentre.com/scriptseqv2)





# LONG OLIGOS UP TO 250 MER

## THE GOLD STANDARD OF OLIGO SYNTHESIS

### SKILLED AND ACCOMPLISHED CAPABILITIES

All Modifications and Oligo Types Synthesized

- Long oligos up to 250 mer
- Fluorescent Molecular Probes
- Ultra-Modified, DNA, RNA, siRNA, Chimeric, Fluorescent, and Antisense Oligos
- Specializing in the design and synthesis of challenging combinations of modifications

Providing all types of oligos for demanding applications and consistent results for more than a decade

Gene Link. Results you can rely on.



toll free: **1-800-GENE LINK** [www.genelink.com](http://www.genelink.com)

190 Saw Mill River Road | Hawthorne, NY 10532  
tel: 914-769-1192 | fax: 914-769-1193








ANNOUNCING

# MOLECULAR CLONING 4



By Michael R. Green, *Howard Hughes Medical Institute, University of Massachusetts Medical School* and Joseph Sambrook, *Peter MacCallum Cancer Institute, Melbourne, Australia*

**M**olecular Cloning: A Laboratory Manual has always been the one indispensable molecular biology laboratory manual for protocols and techniques. The fourth edition of this classic manual preserves the detail and clarity of previous editions as well as the theoretical and historical underpinnings of the techniques presented. Ten original core chapters reflect developments and innovation in standard techniques and introduce new cutting-edge protocols. Twelve entirely new chapters are devoted to the most exciting current research strategies, including epigenetic analysis, RNA interference, genome sequencing, and bioinformatics. This manual is essential for both the inexperienced and the advanced user.

2012, 2,028 pp., illus., appendices, index

Cloth (three-volume set)

\$475

ISBN 978-1-936113-41-5

Paperback (three-volume set)

\$365

ISBN 978-1-936113-42-2

## Contents

### VOLUME 1

#### Part 1 Essentials

1. Isolation and Quantification of DNA
2. Analysis of DNA
3. Cloning and Transformation with Plasmid Vectors
4. Gateway Recombinational Cloning
5. Working with Bacterial Artificial Chromosomes and Other High-Capacity Vectors
6. Extraction, Purification, and Analysis of RNA from Eukaryotic Cells
7. Polymerase Chain Reaction
8. Bioinformatics

### VOLUME 2

#### Part 2 Analysis and Manipulation of DNA and RNA

9. Quantification of DNA and RNA by Real-Time Polymerase Chain Reaction
10. Nucleic Acid Platform Technologies
11. DNA Sequencing
12. Analysis of DNA Methylation in Mammalian Cells
13. Preparation of Labeled DNA, RNA, and Oligonucleotide Probes
14. Methods for In Vitro Mutagenesis

#### Part 3 Introducing Genes into Cells

15. Introducing Genes into Cultured Mammalian Cells
16. Introducing Genes into Mammalian Cells: Viral Vectors

### VOLUME 3

#### Part 4 Gene Expression

17. Analysis of Gene Regulation Using Reporter Systems
18. RNA Interference and Small RNA Analysis
19. Expressing Cloned Genes for Protein Production, Purification, and Analysis

#### Part 5 Interaction Analysis

20. Cross-Linking Technologies for Analysis of Chromatin Structure and Function
21. Mapping of In Vivo RNA-Binding Sites by UV-Cross-Linking Immunoprecipitation (CLIP)
22. Gateway-Compatible Yeast One-Hybrid and Two-Hybrid Assays

#### Appendices

1. Reagents and Buffers
2. Commonly Used Techniques
3. Detection Systems
4. General Safety and Hazardous Material

#### Index



WWW.CSHLPRESS.ORG