KCTCS - UK Pathway to Transfer Guide

Associate in Applied Science – Medical Laboratory Technician

Completion of the following curriculum will satisfy both the requirements for an Associate in Applied Science in MLT at Southeast Community and Technical College as well as UK Core requirements for the University of Kentucky and pre-major requirements for bachelor's degree in Medical Laboratory Sciences.

KCTCS students planning to transfer to the University of Kentucky should work closely with the UK Transfer Advisor each semester to ensure both KCTCS and UK requirements are met.

Summer Session 1	Credit Hrs
MAT 150	3
ENG 101	3
Semester Hours	6

Summer Session 2	Credit Hrs
BIO 135	4
Computer Literacy Course	3
Semester Hours	7

1st Semester		Credit Hrs
MLT 101		3
MLT 112		2
MLT 215		4
PHB 151		2
PHB 152		1
MLT 2781		2.5
BIO 225		4
Semester Hou	urs	18.5

2nd Semester	Credit Hrs
MLT 115	2
MLT 205	3
MLT 216	3
MLT 2782	2.5
CHE 170	4
CHE 175	1
Semester Hours	15.5

Summer Session 3	Credit Hrs
ENG 102	3
STA 220	3
Semester Hours	6

Summer Session 4	Credit Hrs
CHE 180	4
CHE 185	1
Semester Hours	5

3 rd Semester	Credit Hrs
MLT 1191	1.5
MLT 233	3
MLT 225	2
MLT 206	2
MLT 2791	2.5
CHE 270	3
CHE 275	1
Semester Hours	15

4th Semester	Credit Hrs
MLT 1192	1.5
MLT 234	2
MLT 226	2
MLT 2792	2.5
Oral Communications Course	3
BIO 114 or BIO 150	3
BIO 115 or BIO 151	1-2
Semester Hours	15-16

Summer Session 5	Credit Hrs
PSY 110	3
Heritage Course	3
Semester Hours	6

Summer Session 6	Credit Hrs
Social & Behavioral Science Course	3
Humanities	3
Semester Hours	6

Total SCTC Credit Hours: 88-89

KCTCS - UK Pathway to Transfer Guide

Bachelor's in Science – Medical Laboratory Science

5 th Semester	Credit	Hrs
MLS 410	3	
MLS 420	3	
MLS 460	3	
MLS 461	3	
MLS 480	2	
MLS 481	2	
Semester Hours	16	

6 th Semester	Credit Hrs
MLS 430	2
MLS 462	3
MLS 463	3
MLS 464	2
MLS 482	2
MLS 483	2
MLS 471	1
Semester Hours	15

Summer Session 6	Credit Hrs
MLS 440	3
Semester Hours	3

Summer Session 7	Credit Hrs
MLS 450	3
MLS 470	3
Semester Hours	6

Total UK Credit Hours: 40